

# ENGLISH

## A. DEVICE DESCRIPTION

The Ortho-Design VersaTap™, VersaTi™, MiniTi™, MicroTi™ and Duelock™ suture anchors consist of an internal eyelet connected with one, two, or three Force Fibre non-absorbable sutures. The VersaTap™ is a multi-compound anchor made from medical grade PEEK and Titanium, while the VersaTi™, MiniTi™, MicroTi™, and Duelock™ is a single-compound Titanium, Titanium, and PEEK anchor respectively.

## B. INDICATIONS

The VersaTap™ Suture Anchor is intended to be used for soft tissue fixation during general orthopedic surgery. The VersaTap™ Suture Anchor is intended for use in arthroscopic or open surgical approaches for fixation of soft tissue and ligaments to bone/healthy tissue during tendon and ligament repairs, during procedures such as:

- Shoulder: Rotator cuff repair, biceps tenodesis, SLAP repair, Bankart repair
- Foot/Ankle: Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Hallux Valgus
- Elbow: Tennis elbow repair
- Knee: Medial and lateral collateral ligament repair.
- Wrist: Scapholunate ligament reconstruction

The VersaTi™ Suture Anchor is intended for use in arthroscopic or open surgical approaches for fixation of soft tissue and ligaments to bone/healthy tissue during tendon and ligament repairs, during procedures such as:

- Shoulder: Rotator cuff repair, biceps tenodesis, SLAP repair, Bankart repair
- Foot/Ankle: Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Hallux Valgus
- Elbow: Tennis elbow repair
- Knee: Medial and lateral collateral ligament repair, Joint capsule closure
- Wrist: Scapholunate ligament reconstruction
- Hip: Capsular Repair, acetabular labral repair

The MiniTi™/MicroTi™ Suture Anchor is intended to be used for soft tissue fixation during general orthopaedic surgery. The MiniTi™/MicroTi™ Suture Anchor is intended for use in arthroscopic or open surgical approaches for fixation of soft tissue and ligaments to bone/healthy tissue during tendon and ligament repairs, during procedures such as:

- Elbow: Ulnar/Medial Collateral Ligament Repair
- Foot/Ankle: Achilles Tendon Repair, Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Mid-foot Reconstruction, Hallux Valgus Reconstruction, Metatarsal Ligament Repair.
- Hand/Wrist: Scapholunate Ligament Reconstruction and Ulnar/ Radial Collateral Ligament Reconstruction.

The DueLock™ Suture Anchor is intended to be used for soft tissue fixation during general orthopedic surgery. The DueLock™ Suture Anchor is intended for use in arthroscopic or open surgical approaches for fixation of soft tissue and ligaments to bone/healthy tissue during tendon and ligament repairs, during procedures such as:

- Shoulder: Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Capsular Shift and Capsulolabral Reconstruction, Subscapularis Tendon Tears
- Elbow: Ulnar/Radial Collateral Ligament Reconstruction, and Lateral Epicondylitis repair
- Foot/Ankle: Lateral Stabilization, Medial Stabilization, Mid-foot Reconstruction, Hallux Valgus Reconstruction, Metatarsal Ligament Repair, Metatarsal Tendon Repair.

- Hand/Wrist: Scapholunate Ligament Reconstruction and Ulnar/ Radial Collateral Ligament Reconstruction.
- Hip: Acetabular labral repair.

## C. CONTRAINDICATIONS

1. Insufficient quantity or quality of bone.
2. Blood supply limitations and previous infections, which may retard healing.
3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
4. Any active infection or blood supply limitations.
5. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
6. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
7. Do not use for surgeries other than those indicated.

## D. ADVERSE EFFECTS

1. Foreign body reactions.
2. Shoulder dislocation/subluxation
3. Infections (Both deep and superficial)

## E. WARNINGS

1. Caution: Federal law restricts this device to sale by or on the order of a physician.
2. This device is intended to be used by a trained medical professional.
3. An internal fixation device must never be re-used.
4. Do not re-sterilize this device.
5. All metallic implant devices used for this surgical procedure should have the same metallurgical composition.
6. Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight-bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
7. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Ortho-Design delivery system is required for the proper implantation of the device.
8. Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.
9. Detailed instructions on the use and limitations of this device and the brochure ([www.ortho-design.com/recourses](http://www.ortho-design.com/recourses)) should be given to the patient. Your surgeon will guide you in deciding what particular treatment is best for you and explain the benefits, risks, and contraindications associated with the treatment.
10. This is a single-use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
11. Biohazard waste, such as explanted devices, needles, and contaminated surgical equipment, should be safely disposed of in accordance with the institutions' policy.
12. Serious incidents should be reported to Ortho-Design (PTY) Ltd, or an in-country representative, and to the health authority where the incident occurred.

## F. MRI SAFETY INFORMATION

### MR Conditional

*The VersaTap, VersaTi, MicroTi™ and MiniTi™ has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the VersaTap™, VersaTi™, MicroTi™ and MiniTi™ in the MR environment is unknown. Scanning a patient who has this medical device may result in patient injury.*

### MR Safe

#### DueLock™

## G. PRECAUTIONS

1. Surgeons must apply their professional judgment when determining the appropriate suture anchor size based on the specific indication, preferred surgical technique, and patient history.
2. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Ortho-Design provides detailed surgical techniques in print and electronic formats. The Ortho-Design website also provides detailed surgical technique information and demonstrations. Alternatively, contact your Ortho-Design representative for an on-site demonstration.
3. Make sure to use a drill bit or punch to create the bone socket.
4. VersaTap™ suture anchor only: Two light taps with a mallet on the VersaTAP™ insertion handle would identify whether the bone is susceptible to primary self-tapping mode. If the bone is soft enough, make use of the self-tapping (turn 180° forward and 90° back, this allows for bone compression) repeat until the VersaTAP™ is fully seated into the bone.
5. VersaTap™ and Duelock™ suture anchors only: During anchor insertion, ensure that the angle of anchor insertion is coaxial to that of the previously prepared bone socket.
6. VersaTap™ and Duelock™ suture anchors only: Ensure that the anchor body is in full contact with the bone before advancing the anchor body into the prepared bone socket.
7. VersaTi™, MicroTi™ and MiniTi™ suture anchors only: Insertion in very hard bone may require pre-punching a bone socket to avoid damage to the implant.
8. Self-Punching Duelock™ suture anchor only: Ensure that the angle of anchor insertion is perpendicular to the bone.

## H. PACKAGING AND LABELING

1. Ortho-Design devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description, and standard designation number may be found on our website at [www.ortho-design.co.za](http://www.ortho-design.co.za)

## I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions. Certain Ortho-Design instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFO-0023-XX and ANSI/AAMI ST79 for specific information.

## J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials. These devices consist of either one or two components. Each component is manufactured from either titanium alloy or polyetheretherketone (PEEK) or both. **Suture:** See package label for size and type of suture provided with the device. The sutures supplied meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements).

## K. STORAGE CONDITIONS

Non-bioabsorbable devices must be stored in the original unopened packaging, away from moisture, and should not be used after the expiration date.

## L. INFORMATION

Procedures carried out using these devices may be used on the general population. Clinical benefits associated with the use of these devices outweigh the known clinical risks. There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

## VersaTap™, VersaTi™, Duelock™, MiniTi™, MicroTi™, VersaQuad™

MDF-01-06-09

### Suture Anchors:

VersaTap 6mm Suture anchor

VersaTi 5.5mm, 4mm, 3mm Suture anchor

Duelock 2.9mm Knotless/Suture anchor

MiniTi and MicroTi 2mm Suture anchor

VersaQuad - 2 x VersaTap™, PEEK/Ti

Suture Anchors, 6 x 17 mm & 2 x

VersaLat™, PEEK Knotless Anchor

Doc no: ORTHO-

Rev 03



**ORTHO-DESIGN**  
ORTHOPAEDIC INNOVATION

Ortho-Design (PTY) Ltd. Operates from:  
17 Dely Road  
Hazelwood  
Pretoria  
South Africa  
Tel: (+27) 12 807 5172  
[www.ortho-design.com](http://www.ortho-design.com)



TNMC Device Limited  
Krezmina 1,  
10000  
Zagreb  
Croatia  
Tel : +(38) 591 622 9658  
[tnmchr@tnmc.biz](mailto:tnmchr@tnmc.biz)  
[www.tnmc.biz](http://www.tnmc.biz)



This is not a warranty document. For all warranty information, including disclaimers, exclusions, terms, conditions and related provisions, refer to the Ortho-Design website, found at [www.ortho-design.com](http://www.ortho-design.com) whose provisions are incorporated herein by reference.

# ENGLISH

## A. DEVICE DESCRIPTION

The Ortho-Design VersaLat™ 6mm and 7mm knotless anchors consist of a cannulated body with one #2 Force Fiber non-absorbable retention suture passing through. The VersaLat™ is a knotless anchor made from medical grade PEEK, with a Titanium MR-indicating button at the tip.

## B. INDICATIONS

The VersaLat™ Suture Anchor is intended to use for fixation of soft tissue and ligaments to bone/healthy tissue during tendon and ligament repairs, during procedures such as:

- **Shoulder:** Rotator Cuff Repair, Biceps Tenodesis
- **Elbow:** Biceps Tendon Reattachment, Ulnar/Radial Collateral Ligament Reconstruction.
- **Knee:** Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, and Iliotibial Band Tenodesis

## C. CONTRAINDICATIONS

1. Insufficient quantity or quality of bone.
2. Blood supply limitations and previous infections, which may retard healing.
3. Foreign body sensitivity. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
4. Any active infection or blood supply limitations.
5. Conditions that tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
6. The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on skeletally immature patients. The use of this medical device and the placement of hardware or implants must not bridge, disturb or disrupt the growth plate.
7. Do not use for surgeries other than those indicated.

## D. ADVERSE EFFECTS

1. Foreign body reactions.
2. Shoulder dislocation/subluxation
3. Infections (Both deep and superficial)

## E. WARNINGS

1. Caution: Federal law restricts this device to sale by or on the order of a physician.
2. This device is intended to be used by a trained medical professional.
3. An internal fixation device must never be re-used.
4. Do not re-sterilize this device.
5. Postoperatively and until healing is complete, fixation provided by this device should be considered as temporary and may not withstand weight-bearing or other unsupported stress. The fixation provided by this device should be protected. The postoperative regimen prescribed by the physician should be strictly followed to avoid adverse stresses applied to the device.
6. Pre-operative and operating procedures, including knowledge of surgical techniques and proper selection and placement of the device, are important considerations in the successful utilization of this device. The appropriate Ortho-Design delivery system is required for the proper implantation of the device.
7. Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Device removal should be followed by adequate postoperative management.

8. Detailed instructions on the use and limitations of this device and the brochure ([www.ortho-design.com/recourses](http://www.ortho-design.com/recourses)) should be given to the patient. Your surgeon will guide you in deciding what particular treatment is best for you and explain the benefits, risks, and contraindications associated with the treatment.
9. This is a single-use device. Reuse of this device could result in failure of the device to perform as intended and could cause harm to the patient and/or user.
10. Biohazard waste, such as explanted devices, needles, and contaminated surgical equipment, should be safely disposed of in accordance with the institutions' policy.
11. Serious incidents should be reported to Ortho-Design (PTY) Ltd, or an in-country representative, and to the health authority where the incident occurred.

## F. MRI SAFETY INFORMATION

### MR Conditional

*Not evaluated.*

### MR Safe

*Not evaluated.*

## G. PRECAUTIONS

1. Surgeons must apply their professional judgment when determining the appropriate suture anchor size based on the specific indication, preferred surgical technique, and patient history.
2. Surgeons are advised to review the product-specific surgical technique prior to performing any surgery. Ortho-Design provides detailed surgical techniques in print and electronic formats. The Ortho-Design website also provides detailed surgical technique information and demonstrations. Alternatively, contact your Ortho-Design representative for an on-site demonstration.
3. Make sure to use a drill bit or punch to create the bone socket.
4. Make sure to use a tap/drill for both the 6mm and 7mm VersaLat. This is essential due to the difference in size and indication purposes.
5. During anchor insertion, ensure that the angle of anchor insertion is coaxial to that of the previously prepared bone socket.
6. Ensure that the anchor body is in full contact with the bone before advancing the anchor body into the prepared bone socket.

## H. PACKAGING AND LABELING

1. Ortho-Design devices should be accepted only if the factory packaging and labeling arrive intact.
2. Contact Customer Service if the package has been opened or altered.
3. All of the symbols used on the labeling along with the title, description, and standard designation number may be found on our website at [www.ortho-design.co.za](http://www.ortho-design.co.za)

## I. STERILIZATION

This device is provided sterile. Check the package labeling for more information. This device should never be re-sterilized under any conditions. Certain Ortho-Design instruments that may be used during this procedure are provided non-sterile and must be adequately cleaned and sterilized prior to use or re-use. Please refer to DFU-0023-XX and ANSI/AAMI ST79 for specific information.

## J. MATERIAL SPECIFICATIONS

Refer to the package label for the materials. These devices consist of either one or two components. Each component is manufactured from either titanium alloy or polyetheretherketone (PEEK) or both. **Suture:** See package label for size and type of suture provided with the device. The sutures supplied meet or exceed U.S. and European Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements).

## K. STORAGE CONDITIONS

Non-bioabsorbable devices must be stored in the original unopened packaging, away from moisture, and should not be used after the expiration date.

## L. INFORMATION

- Procedures carried out using these devices may be used on the general population.
- Clinical benefits associated with the use of these devices outweigh the known clinical risks.
- There are no unacceptable residual risks or uncertainties associated with the clinical use of these devices.

**VersaLat™**

Doc no: ORTHO-MDF-01-06-10

**Knotless Anchor:**

*VersaLat 7.5mm, 6.5mm, 5.5mm Knotless anchor*

Rev 03



**ORTHO-DESIGN**  
ORTHOPAEDIC INNOVATION

Ortho-Design (PTY) Ltd. Operates from:  
17 Dely Road  
Hazelwood  
Pretoria  
South Africa  
Tel: (+27) 12 807 5172  
[www.ortho-design.com](http://www.ortho-design.com)



TNMC Device Limited  
Krezmina 1,  
10000  
Zagreb  
Croatia  
Tel : +(38) 591 622 9658  
[tnmchr@tnmc.biz](mailto:tnmchr@tnmc.biz)  
[www.tnmc.biz](http://www.tnmc.biz)



This is not a warranty document. For all warranty information, including disclaimers, exclusions, terms, conditions and related provisions, refer to the Ortho-Design website, found at [www.ortho-design.com](http://www.ortho-design.com) whose provisions are incorporated herein by reference.