




## Duelock™

Knotted Suture Anchor  
Engineered for optimal  
pull-out strength



### Shoulder

#### 2.9mm PEEK Duelock™



- one #2 UHMWPE suture
- single-use inserter
- bankart repair, labral repair, SLAP repair, elbow surgery



### Hip


#### 2.9mm PEEK Duelock™

- one #2 UHMWPE suture
- single-use inserter
- hip labral repair




## Features and Benefits

- Provides the surgeon with a knotted fixation method
- Increased pull-out strength and a reduced footprint
- A tapered thread design to maximize cortical compression



## Materials & Specifications

- Anchor Material: 100% Medical grade Polyether ether ketone (PEEK)
- Anchor Size: 2.9mm x 10mm
- Suture Material: One #2 Force Fibre UHMWPE braided suture (non-absorbable)
- Introducer: Single-use, cannulated introducer composing of medical grade stainless steel and silicone



## Extensive Application Range

- The DueLock™ ancillary can be adjusted for shoulder, elbow and hip surgery

