




VersaLat™

*Outstanding suture retention
Ultimate Knotless fixation
Increased functionality*



6mm PEEK VersaLat™


- one #2 UHMWPE suture
- single-use inserter

7mm PEEK VersaLat™

- one #2 UHMWPE suture
- single-use inserter

8mm PEEK VersaLat™

- one #2 UHMWPE suture
- single-use inserter




T: 012 807 5172
www.ortho-design.com
info@ortho-design.co.za

ORTHO-DESIGN
ORTHOPAEDIC INNOVATION

Features and Benefits

- A one-handed technique promotes ease-of-use and decreases operation time
- Allows knotless fixation of a large number of sutures due to the adjustable suture loop
- Unique locking mechanism for effective & efficient suture tensioning
- Eliminates knots associated with tissue irritation
- A retention suture that allows for the fixation of any remaining ligament fragments
- Engineered counter thread allows for independent rotation of the inner shaft upon anchor insertion



Extensive Application Range

- The 6mm VersaLat™ to be used as a knotless anchor in single and double-row rotator cuff repair surgery
- The 7mm and 8mm VersaLat™ to be used for biceps tenodesis repair
- Can be used as a general interference screw or internal brace

Materials & Specifications

- Anchor Material: 100% Medical grade Polyether Ether Ketone (PEEK)
- Retention Suture Material: One #2 Force Fibre UHMWPE braided suture (non-absorbable)
- Introducer: Single-use, cannulated introducer composing of medical grade stainless steel and silicone

