

Free Needles Suture Passing



FREE NEEDLES

Free needles for manipulation of sutures during soft tissue fixation surgery.

FEATURES AND BENEFITS

- Allows seamless manipulation of sutures through soft tissue without the need for pre-attached needles
- Can be used for all suture types and sizes



MATERIALS & SPECIFICATIONS

NMO-06, 26mm Taper Point Needle

USP 3/0 Nylon Suture

18mm (9mm Loop)

EXTENSIVE APPLICATION RANGE

Can be used for all suture types and sizes

References

FN13S	Free needle sterile small 13mm
FN23M	Free needle sterile medium 23mm
FN40L	Free needle sterile large 40mm