


ADJUSTABLE CORTICAL FIXATION LOOP



The Adjustable Cortical Fixation Loop device offers leading strength and incorporating a unique cradle design to help protect the graft during loop reduction.

Ideal for primary or auxiliary fixation ACL or PCL reconstruction techniques. Preloaded with UHMWPE suture and flipping suture for added procedure efficiency.

FEATURES AND BENEFITS

- **High tensile strength.**
- **Graft support frame with large contact area.**
- **Unique cradle design to help protect the graft during loop reduction.**
- **Enhanced precision.**
- **Quick and secure fixation.**

MATERIALS & SPECIFICATIONS

Loop: 15mm - 60mm


Sutures: UHMWPE Sutures

Plate: Medical Grade Titanium (Ti-6Al-4V)

EXTENSIVE APPLICATION RANGE

Use in the fixation of bone and soft tissue in orthopedic procedures requiring ligament or tendon reconstruction.

Can be used in Single-Bundle soft tissue fixation; Double-Bundle soft tissue fixation; Bone-Tendon-Bone fixation.



References

OCBGEFRB1U Adjustable Cortical Fixation Loop Device
OCBGEPLB1U Fixation Loop Extended Plate