


VersaPEEK™

SUTURE ANCHOR



The VersaPEEK™ is a 100% PEEK suture anchor used in rotator cuff, general soft tissue fixation, and complex fracture surgeries. The anchor accommodates a wide variety of suture or suture tape configurations and is engineered for ultimate pull-out strength.

ADVANTAGES:

- Minimal post-operative imaging effects
- Engineered for ultimate mechanical properties (pull-out strength, tensile strength, etc.)
- Variety of suture/suture tape configurations
- Each size comes in a closed and open-handle design - Needles are retained inside the handle.

INDICATIONS:

- Rotator Cuff Repair.
- Foot and Ankle Ligament Stabilization and Tendon Repairs.
- Ligament and Tendon Repair in the Elbow.
- Biceps Repairs and Tenodesis.

COMPOSITION:

Anchor body:
100% PEEK

Suture Material:
UHMWPE braided, non-absorbable

Needles:
MO-6 26mm half-round taper point



References

APP001	VersaPEEK, PEEK Suture Anchor, 5.5 x 15mm
APP001-N	VersaPEEK, PEEK Suture Anchor, 4.5 x 15mm, w/needles
APP002	VersaPEEK, PEEK Suture Anchor, 4.5 x 15mm
APP002-N	VersaPEEK, PEEK Suture Anchor, 4.5 x 15mm, w/needles